MicroRNA

Subjects: Biotechnology & Applied Microbiology

Contributor: Paola Campomenosi , Francesca Precazzini , Michela Alessandra Denti

MicroRNA (miRNAs) are short, single-stranded ncRNAs that can be measured in biofluids as minimally invasive biomarkers of disease, including cancer. Their reliable measurement however is challenging, therefore several methods are being developed for miRNA quantitation. Here the technologies currently in use or being developed for quantitation of miRNAs are described.

MicroRNA biofluids methods biomarkers cancer

1. Introduction

Blood circulating microRNAs (c-miRNAs) are considered valuable biomarkers since they have show high stability, resistance to ribonucleases and to severe physicochemical conditions in body fluids. This favours their use in clinical applications. Differential expression of c-miRNAs showed big potential for cancer screening. Thus, these small RNA molecules may function as clinical biomarkers to identify the presence of tumors, to select treatment strategies and to predict outcomes as explained in Figure 1. Composition and levels of extracellular miRNAs can vary as a consequence of diseases or other conditions ^[1]. Several reports show that miRNAs are sensitive to physiological or pathological changes, specific to the pathology of interest, and they represent a reliable indicator of the presence of disease before clinical symptoms appear ^[2]. Moreover miRNAs have shown great potential as diagnostic biomarkers for neurodegenerative disorders, cardiovascular diseases, endometriosis, cancer and several other diseases.



Figure 1. miRNAs are involved in all cellular processes and dysregulated in disease, potentially representing optimal biomarkers. They can be purified from different tissues or liquid biopsies and then measured with the different technologies described in the text. Differential expression of miRNAs may indicate an increased risk of disease or have a diagnostic or prognostic value, leading to clinical decisions.

2. Technical Limitations in Measurement of c-miRNAs as Biomarkers

If biological complexity is already difficult to capture, there are also technical limitations in c-miRNAs measurement that make it difficult to develop a diagnostic or prognostic test directly translatable in clinical settings. In fact, despite the growing interest of the scientific community, proved by the enormous amount of articles published in recent years (a search for "circulating microRNAs AND cancer biomarker" on PubMed, performed in September 2020, yielded 123 entries in 2012, 189 in 2013, 254 in 2014, 335 in 2015, 329 in 2016, 338 in 2017, 356 in 2018 and 310 in 2019) only two miRNA based diagnostic tests (Biomild, Cosmos II) have reached phase 4 in clinical trials ^[3]. The main reasons for the limited clinical application may be classified into three categories: 1) pre-analytical, 2) analytical and 3) post-analytical factors. These factors affect measurement accuracy and reproducibility independently of its clinical application.


Table 1. Summary of the different technical limitations during miRNA quantification.

Phase of the Analysis	Challenges					
	 Gender, smoking status and physiological conditions (diet, active exercise, circadian rhythm) 					
Pre-analytical phase	Comorbidities					
	[–] Type of sample (plasma, serum, saliva, urine, or FFPE and flash-frozen tissues)					
	 Modality of sample collection (e.g., venipuncture) or time from collection to processing 					
	 For plasma: Which type of anticoagulant 					
	 Sample processing (centrifugations steps) 					
	 Sample storage 					
Analytical phase	RNA extraction-different kits					
	[–] Primer design: Difficulties to discriminate pri-miRNA, pre-miRNA and mature forms					

Phase of the Analysis	Challenges					
	 Reverse transcription: Specific primers or universal primers could introduce amplification biases 					
	Extraction/quantification efficiency affected by miRNA sequence composition					
Post-analytical phase (normalization)	 Single reference molecule 					
	 Set of reference molecules 					
	 Mathematical normalization 					
	 Different software for the selection of the best reference molecule (Normfinder, geNorm, Bestkeeper) 					
	 Calibration curve 					

3. Indirect Methods for miRNA Measurement, Relying on RNA Extraction & Reverse Transcription (RT)

Several methods have been applied to quantify c-miRNAs, including quantitative real-time PCR, digital PCR, microarray, and high-throughput small RNA-sequencing. These methods need RNA isolation, in some instances small RNA enrichment, and/or reverse transcription before proceeding to miRNA quantification. A summary of the different techniques described in this review is depicted in Figure 2.





4. Direct Methods for miRNA Measurement from Biofluids

HTG EdgeSeq

In 2018 an approach called HTG EdgeSeq, developed by HTG Molecular Diagnostic, Inc. (Tucson, AZ, USA), was described for the direct screening of plasma c-miRNAs. By eliminating RNA extraction and reverse transcription, it

simplified both library and sample preparation for RNA and miRNA targeted sequencing. Probes are used to form double stranded molecules with targets and protect the complexes from degradation with S1 nuclease. The probes also carry "wings", extensions that will allow introduction of adapters for NGS by PCR, protected by "wingmen". Single stranded nucleic acids (targets and unbound probes), will be degraded, leaving the target-probe complexes available for subsequent steps (Figure 3). This approach allows direct quantification of transcripts or miRNAs, allowing to increase the throughput and to considerably reduce bioinformatic analyses and time to results (~36 hours). Profiles of more than 2000 miRNA transcripts can be performed starting from a relatively small amount of initial sample. Normalization can be performed as classical sRNA-Seq. Moreover, the HTG EdgeSeg costs are consistent with other next-generation sequencing (NGS) approaches. In a study by Songia and coworkers peripheral blood from four healthy subjects was collected. Samples were barcoded, purified and quantified using a KAPA Library Quantification kit 4. Overall HTG-sequencing with Illumina NextSeq was performed leading to the evaluation of the expression of 2083 miRNAs starting from just 15 µL of plasma. In addition, the study compared the performances of techniques that are commonly used for miRNA detection like guantitative PCR (gPCR) and the chip-based digital PCR (dPCR) and a strong correlation between gPCR and dPCR with HTG was found ^[4]. In another comparative study, the performance of HTG Edgeseq was compared with that of classical RNA-seq, Fireplex and Nanostring nCounter platforms in measurement of synthetic miRNA pools or miRNAs from plasma samples [146]. HTG Edgeseg showed the smallest coefficient of variation when repeated measures were performed and low bias in measuring "ratiometric pools" of synthetic miRNAs. Moreover, it could detect variation in the levels of placenta-specific miRNAs between pregnant and non-pregnant women ^[5].




Figure 3. Diagram of HTG Edgeseq workflow. The hybridization, S1 nuclease treatment and PCR steps are depicted. These steps allow it to protect probes bound to target miRNAs, removing other RNAs and unbound probes. Wash steps are not shown.

Table 2. Summary of the different time, throughput and multiplexing of the different methods described in this entry.

	Method	Time	ThroughputMultiplexing		
INDIRECT MEASUREMENTS	RNA extraction and RT	qPCR	<6 h	Medium	Low
		digital PCR	<6 h	Medium	Low
		Microarray	~2 days	Low	High
		NGS	1–2 weeks	Medium	High

	Method	Time	Throughput		
	RT only	RT-qPCR-DS	<4 h	Medium	Low
		Cell to Ct RT- qPCR			
		Direct S-Poly(T) Plus	<4 h	Medium	Low
		RT-qPCR	<4 h	Medium	Low
	RNA extraction only	Nanostring Platform	~2 days	Low	High
DIRECT MEASUREMENTS		HTG EdgeSeq	~36 h	Medium	High
		bDNA	~12 h	Low	High
		Fireplex	~6 h	High	Medium
		Chem-NAT	~3 h	Medium	Low

5. Conclusions

A comprehensive knowledge of c-miRNAs and their origin, as well as of other classes of non-coding RNAs including IncRNAs and circular-RNAs, has still to be obtained. The combination of miRNAs with other biomarkers should be evaluated. However, a clinical diagnostic test requires not only specific analytes that can identify the presence of a condition, be it cancer or other diseases, but also a reliable method for their quantification. The present review describes up to date techniques that can be applied, and whose performance can be compared, to address the issue of reliable miRNA biomarker quantification. With stringent guidelines, a detailed description of the procedures used for miRNA quantification in published papers and increased application of direct measurement methods, we expect that a step closer to their clinical application can be reached. Independent validation studies are required before microRNA measurement can be translated into clinical applications.

References

- 1. Li, Y.; Kowdley, K.V. MicroRNAs in common human diseases. Genom. Proteom. Bioinform. 2012, 10, 246–253.
- De Guire, V.; Robitaille, R.; Tétreault, N.; Guérin, R.; Ménard, C.; Bambace, N.; Sapieha, P. Circulating miRNAs as sensitive and specific biomarkers for the diagnosis and monitoring of human diseases: Promises and challenges. Clin. Biochem. 2013, 46, 846–860.

- Seijo, L.M.; Peled, N.; Ajona, D.; Boeri, M.; Field, J.K.; Sozzi, G.; Pio, R.; Zulueta, J.J.; Spira, A.; Massion, P.P.; et al. Biomarkers in Lung Cancer Screening: Achievements, Promises, and Challenges. J. Thorac. Oncol. 2019, 14, 343–357.
- 4. Songia, P.; Chiesa, M.; Valerio, V.; Moschetta, D.; Myasoedova, V.A.; D'Alessandra, Y.; Poggio, P. Direct screening of plasma circulating microRNAs. RNA Biol. 2018, 15, 1268–1272.
- Godoy, P.M.; Barczak, A.J.; DeHoff, P.; Srinivasan, S.; Etheridge, A.; Galas, D.; Das, S.; Erle, D.J.; Laurent, L.C. Comparison of Reproducibility, Accuracy, Sensitivity, and Specificity of miRNA Quantification Platforms. Cell Rep. 2019, 29, 4212–4222.e5.

Retrieved from https://encyclopedia.pub/entry/history/show/18361