Characterization of the Oxidation Products of Lipids


Subjects: Geochemistry & Geophysics Contributor: Rontani Jean-François

Phototrophic organisms (mainly green plants, algae, cyanobacteria and some protists) carry out photosynthesis that is, conversion of sunlight energy, carbon dioxide and water into organic materials. Due to the generation of highly reactive oxygen species (ROS) during photosynthetic electron transport, these organisms are particularly sensitive to oxidative damages. Lipids (hydrocarbons, pigments, terpenoids, free fatty acids, acylglycerides, phospholipids, galactolipids, cutins, suberins and waxes) are important components of phototrophic organisms, accounting for 16–26% of organic content in phytoplankton and up to 45% in the green alga *Botryococcus Braunii*. The applications of gas chromatography-mass spectrometry techniques for the characterization of photooxidation and autoxidation products of unsaturated lipids of senescent phototrophic organisms will be introduced.

Keywords: senescent phototrophs ; unsaturated lipids ; Tracers ; abiotic oxidation ; GC-MS


1. Chlorophyll Phytyl Side-Chain

Attack of ¹O2 on the tri-substituted double bond of the chlorophyll phytyl side-chain affords two allylic hydroperoxides, which may be recovered in the form of 6,10,14-trimethylpentadecan-2-ol and 3-methylidene-7,11,15-trimethylhexadecan-1,2-diol (phytyldiol) after NaBH4 reduction and alkaline hydrolysis ^[1] (Scheme 1). The stable and highly specific phytyldiol was proposed as biogeochemical marker of chlorophyll photodegradation in the natural environment ^[2]. In contrast, free radical oxidation (autoxidation) of chlorophyll phytyl side-chain and subsequent reduction and hydrolysis gives 3,7,11,15-tetramethyl-hexadec-2-en(Z/E)-1,4-diols and 3,7,11,15-tetramethyl-hexadec-1-en-3-ol (isophytol) ^{[3][4]} (Scheme 1). These compounds have been proposed as specific tracers of chlorophyll phytyl side-chain autoxidation in environmental samples ^{[3][4]}.



Scheme 1. Photooxidation and autoxidation of chlorophyll phytyl side-chain.

TOF mass spectra of the TMS derivatives of phytyldiol and 3,7,11,15-tetramethylhexadec-3-en(Z/E)-1,2-diols show intense and specific fragment ions at m/z 353.3235 resulting from classical α -cleavage between the carbon atoms 1 and 2 bearing the two TMS ether groups ^[5], while the spectra of TMS derivatives of 3,7,11,15-tetramethyl-hexadec-2-en(Z/E)-1,4-diols are dominated by a fragment ion at m/z 245.1388 corresponding to α -cleavage between carbon atoms 4 and 5 (Scheme 2).



Scheme 2. Main EI mass fragmentations of TMS derivatives of phytyldiol, 3,7,11,15-tetramethylhexadec-3-en(Z/E)-1,2-diols and 3,7,11,15-tetramethyl-hexadec-2-en(Z/E)-1,4-diols.

Note that the loss of a methyl radical by the molecular ion of TMS derivatives of phytol and isophytol also affords a fragment ion at m/z 353.3235. Monitoring ions at m/z 353.3235 and 245.1388 thus allows simultaneous characterization and quantification of phytol and its main photooxidation and autoxidation products in natural samples (see example given in **Figure 1**).





Figure 1. Partial TOF ion chromatograms (m/z 353.3235 and 245.1388) showing the presence of TMS derivatives of phytol and its main photooxidation and autoxidation products in senescent cells of the diatom Thalassiosira sp.

2. Δ⁵-Sterols


Reaction of ¹O2 with the double bond of Δ^5 -sterols mainly affords a Δ^6 -5 α -hydroperoxide and to a lesser extent Δ^4 -6 α/β -hydroperoxides ^{[6][7]} (Scheme 3). Under environmental conditions Δ^6 -5 α -hydroperoxide undergoes fast allylic rearrangement to unstable and unspecific 7 α/β -hydroperoxides (Scheme 3). Δ^4 -Stera-3 β ,6 α/β -diols resulting from NaBH4-reduction and alkaline hydrolysis of Δ^4 -6 α/β -hydroperoxides were thus proposed as specific tracers of type II photosensitized oxidation of the corresponding Δ^5 -sterols ^{[8][9]}.



Scheme 3. Photooxidation and autoxidation of Δ^5 -sterols.

Autoxidation of Δ^5 -sterols mainly affords unstable and unspecific $7\alpha/\beta$ -hydroperoxides after hydrogen atom abstraction at the allylic carbon atom 7 ^[10]) (<u>Scheme 3</u>). Smaller proportions of isomeric $5\alpha,6\alpha$ - and $5\beta,6\beta$ -epoxysterols are also produced after addition of peroxyl radical to the double bond ^[11] (<u>Scheme 3</u>). Stable and specific $3\beta,5\alpha,6\beta$ -trihydroxysterols resulting from the hydrolysis of these epoxides during alkaline hydrolysis and in environmental conditions were proposed as specific tracers of the autoxidation of Δ^5 -sterols ^{[8][9]}.

TOF mass spectra of Δ^4 -stera-3 β ,6 α / β -diol TMS derivative exhibit an intense and interesting fragment ions at [M – 143.0887]⁺ resulting from double bond ionization and subsequent hydrogen migrations and cleavages of the C1–C10 and C4–C5 bonds ^[12] (Scheme 4). Due to steric hindrance, the classical silvlation reagents only silvlate 3 β ,5 α ,6 β -trihydroxysterols to their 3 and 6 positions ^[13], and during ionization their TMS derivatives very easily lose a neutral molecule of water and thus exhibit mass spectra that are very similar to those of Δ^4 -stera-3 β ,6 β -diol TMS derivatives (Scheme 4).



Scheme 4. Proposed formation pathways of the fragment ion $[M - 143.0887]^+$ in TOF mass spectra of Δ^4 -stera-3 β ,6 β -diol and 3 β ,5 α ,6 β -trihydroxysterol TMS derivatives.

3. Unsaturated Fatty Acids

Type II photosensitized oxidation and autoxidation rates of unsaturated fatty acids logically increase with their number of double bonds $^{[14][15]}$. Unfortunately, oxidation products of the more reactive polyunsaturated fatty acids (PUFAs) are not sufficiently stable under environmental conditions to be used as tracers of these degradation processes in situ. Note that isoprostanoids (cyclopentane-containing oxylipins) resulting from autoxidation of C18, C20 and C22 PUFAs are often used as biomarkers for in vivo oxidative stress in animals and plants $^{[16]}$. These compounds could be detected in higher plants and algae by using GC–MS in negative-ion chemical ionization (NICI) mode $^{[16][17]}$. However, in the literature there are no reports of such compounds in environmental samples.

GC-QTOF allows a clean characterization and quantification of TMS derivatives of MUFA oxidation products in autotrophic organisms and environmental samples. **Figure 2** gives some examples of the technique application showing typical profiles of visible light-induced, (visible + UV) light-induced and autoxidative degradation products.




Figure 2. Partial TOF ion chromatograms showing TMS derivatives of oleic acid oxidation products in senescent cells of the haptophyte Emiliania huxleyi irradiated by visible light (**A**), and after aging (**B**), and of palmitoleic acid in the diatom Thalassiosira sp. irradiated by (visible + UV) light (**C**).

MRM analyses of TMS derivatives of MUFA oxidation products involve intense and selective transitions from the ions resulting from α -cleavage (precursor ions) to the fragment ion at m/z 129 (product ion). Note that this transition is more efficient with precursor ions containing the terminal methyl group than with precursor ions containing the TMS ester group, which can easily lose neutral TMSOH molecules (**Figure 3**). **Figure 4** gives an example of how MRM analyses can be applied.





Figure 3. Low energy collision-induced dissociation (CID)-MS/MS (5 eV) of fragment ions at m/z 227 and 329.



Retention time (min)


Figure 4. MRM chromatogram (m/z 199 \rightarrow 129, m/z 213 \rightarrow 129, m/z 329 \rightarrow 149 and m/z 343 \rightarrow 163) showing the presence of TMS derivatives of palmitoleic acid (C16:1 ω 9) oxidation products in senescent cells of Thalassiosira sp. irradiated by sunlight.

4. Pentacyclic Triterpenes

Pentacyclic triterpenes and their derivatives, which are widely found in angiosperms ^[18], are divided into three main classes, that is, lupanes, oleananes and ursanes.


4.1. Lupanes

Type-II photooxidation and autoxidation of lupanes have so far only been studied for betulin $\frac{[19]}{1}$, but the results obtained can be extended to lupeol or betulinic acid (the main triterpenoids with betulin of the lupane group). ¹O2 reacts slowly with the C20–C29 double bond of betulin and specifically produces lup-20(30)-ene-3 β ,28-diol-29-hydroperoxide, which can be quantified after NaBH4 reduction in the form of lup-20(30)-ene-3 β ,28,29-triol (<u>Scheme 5</u>). Lup-20(30)-ene-3 β ,28,29-triol, lup-20(30)-ene-3 β ,29-diol (arising from lupeol) and lup-20(30)-ene-3 β ,29-diol-28-oic acid (arising from betulinic acid) constitute useful specific tracers of photooxidation of lupanes in angiosperms.



Scheme 5. Photooxidation and autoxidation of lupanes.

The EI mass spectra of the TMS derivatives of lup-20(30)-ene-3 β ,28,29-triol and lupan-20-one-3 β ,28-diol exhibit intense fragment ions at m/z 481.3860 and m/z 395.3305, respectively, whose formation involves elimination of a neutral molecule of TMSOH and subsequent loss of the CH2OTMS group borne by the carbon 28 ^[19] (Scheme 6). These fragment ions make good candidates for monitoring type-II photosensitized oxidation and autoxidation of betulin, respectively, in environmental samples. As the formation of these ions involves the loss of the group borne by carbon 28, they can be also used as tracers of the oxidation of lupeol and betulinic acid. **Figure 5** gives an example of the specific fragment ion at m/z 395.3305 applied for monitoring lupane autoxidation in environmental samples.



Scheme 6. Main EI fragmentations of the TMS derivatives of lup-20(30)-ene-3β,28,29-triol and lupan-20-one-3β,28-diol.





Figure 5. Partial TOF ion chromatograms (m/z 395.3305, 498.3908 and 500.4053) showing the presence of autoxidation products of lupeol and betulin in higher plant debris collected in the Rhône River.

4.2. Ursanes and Oleanes

Studies on type II photosensitized oxidation and autoxidation of ursanes and oleanes have mainly focused on α - and β amyrins ^[20]. α - and β -amyrins were found to be totally unaffected during photodegradation experiments, due to steric hindrance preventing ¹O2 reaction with their double bond ^[20]. Autoxidation of amyrins mainly involves hydrogen abstraction and specifically produces 11 α -hydroperoxyamyrins ^[20] (Scheme 7).



Scheme 7. Autoxidation of α - and β -amyrins.

These hydroperoxides, which appeared to be unaffected by NaBH4 reduction, are thermally cleaved to the corresponding 11-oxoamyrins during GC or GC-MS analyses using hot injectors (<u>Scheme 7</u>). 11-Oxoamyrins are sufficiently stable and specific to serve as tracers of amyrin autoxidation in senescent angiosperms or environmental samples.

Specific fragment ions at m/z 512.4063 [M]^{+*}, 383.3308, 273.2213 and 232.1822 appeared to be useful for GC-QTOF monitoring of TMS derivatives of 11-oxo-amyrins in environmental samples (see example given in **Figure 6**A). MRM analyses using the transitions m/z 273 \rightarrow 135 and m/z 232 \rightarrow 217 also appeared to be well suited to the detection of traces of these compounds (see **Figure 6**B).




Figure 6. Partial TOF ion chromatogram (m/z 512.4063, 383.3308, 273.2213 and 232.1822) (**A**) and MRM chromatogram (m/z 273 \rightarrow 135 and m/z 232 \rightarrow 217) (**B**) showing the presence of oxidation products of amyrins in senescent leaves of Quercus ilex.

5. Dehydroabietic Acid

Fragment ions at m/z 191.0887, 234.1435 and 237.1638 can be used in GC-QTOF analyses to characterize TMS derivatives of the main autoxidation products of dehydroabietic acid ($7\alpha/\beta$ -hydroxydehydroabietic acids). However, MRM analyses using the highly specific transitions m/z 234 \rightarrow 191 and m/z 460 \rightarrow 417 [M – isopropyl group]⁺ emerged as better suited to detecting traces of these compounds in environmental samples (**Figure 7**).





Figure 7. Partial MRM chromatogram (m/z $234 \rightarrow 191$ and m/z $252 \rightarrow 237$) showing the presence of autoxidation products of dehydroabietic acid in senescent needles of Pinus halepensis.


6. Highly Branched Isoprenoid (HBI) Alkenes

HBI alkenes (exhibiting 1-6 double bonds) are produced by some marine and freshwater diatoms belonging to the Berkeleya, Haslea, Navicula, Pleurosigma, Pseudosolenia and Rhizosolenia genera [21][22]. During the senescence of these organisms, ¹O2 attack is focused on the lesser sterically-hindered trisubstituted double bonds of these alkenes affording 2 or 4 allylic hydroperoxides according to the E or Z configuration of the double bond ^[23]. As an example, Scheme 8 shows type II photosensitized oxidation of Z and E isomers of HBI III, which are ubiquitous throughout the world's oceans ^[24]. In this case, ¹O2 attack acts mainly on the C9–C10 double bond and to a lesser extent to the more sterically-hindered C7-C20 double bond affording 9- and 7-hydroperoxides, respectively, as the major oxidation products. Autoxidation processes also act very quickly on HBI III, producing numerous autoxidation products, but predominantly 9-hydroperoxides resulting from hydrogen atom abstraction at the allylic carbon 11 (Scheme 8) [25]. Indeed, the major oxidation pathway of this compound involves hydrogen abstraction at the bis-allylic C8 position to afford conjugated dienes, which are particularly prone to peroxyl radical additions and readily undergo copolymerization with oxygen (Scheme 8). Consequently, the 7-alcohol resulting from NaBH4-reduction of the corresponding hydroperoxide could be used as specific tracer of type II photosensitized oxidation of HBI III (Scheme 8). However, the reduction products of 9-hydroperoxides will only be indicative of oxidation of this specific HBI alkene. Unfortunately, in the case of HBI alkenes (such as HBI III) possessing several trisubstituted double bonds, photooxidation and autoxidation products are unable to accumulate due to the involvement of fast secondary oxidation reactions [25]. All these tracers can thus only serve to give qualitative indications.



Scheme 8. Type II photosensitized oxidation and autoxidation of HBI III.

EI mass spectra of the TMS derivatives of the 9-alcohols resulting from HBI III oxidation exhibit an intense fragment ion at m/z 213.1670 corresponding to α -cleavage relative to the TMS ether group ^[26] (Scheme 9). This fragment ion can readily lose a neutral molecule of TMSOH to give a fragment ion at m/z 123.1170, or undergo a hydrogen transfer with concerted cleavage of the bond between carbon atoms 3 and 4, yielding a fragment ion at m/z 143.0887 (Scheme 9). In the case of the 7-alcohol, α -cleavage relative to the TMS ether group affords two fragment ions at m/z 295.2452 and 321.2610, which are then cleaved in the α position relative to the ionized TMS ether group after hydrogen transfers to give fragment ions at m/z 183.1201 and 181.1044, respectively (Scheme 9).



Scheme 9. Main EI mass fragmentations of TMS derivatives of 9- and 7-alcohols resulting from oxidation of HBI III.

Oxidation products of HBI III were only characterizable in environmental samples in MRM mode using the m/z 213 \rightarrow 123, m/z 213 \rightarrow 143, m/z 295 \rightarrow 183 and m/z 321 \rightarrow 181 transitions ^{[25][26][27]}. An applied example is given in **Figure 8**.





Figure 8. Partial MRM chromatograms (m/z 213 \rightarrow 123, m/z 213 \rightarrow 143, m/z 321 \rightarrow 181 and m/z 295 \rightarrow 183) showing the presence of oxidation products of HBI III in diatoms collected in Commonwealth Bay (Antarctic).

7. Alkenones

Alkenones are a class of mono-, di-, tri-, tetra- and penta-unsaturated C35–C40 methyl and ethyl ketones, which are produced by certain haptophytes ^{[28][29][30][31][32]}. The unsaturation ratio of C37 alkenones, which is defined by the equation: UK'37 = [C37:2]/([C37:2] + [C37:3]) (where [C37:2] and [C37:3] are the concentrations of di- and tri-unsaturated C37 methyl alkenones, respectively) varies positively with the growth temperature of the alga ^{[33][34]} and is thus now routinely used for paleotemperature reconstructions (e.g., ^{[35][36]}). Due to the trans- geometry of the alkenone double bonds ^[37], which is poorly reactive with ¹O2 ^[38]), alkenones are not affected by type II photosensitized oxidation processes ^{[39][40]}. However, they are highly reactive to autoxidation processes ^[41]. Autoxidation of alkenone double bonds (separated by five or three methylene groups) affords six hydroperoxides as in the case of MUFAs . Isomeric alkenediols resulting from NaBH4 reduction of these oxidation products could make very useful indicators of autoxidative alterations of the unsaturation ratio UK'37, but unfortunately they fail to accumulate due to the subsequent oxidation of the other double bonds ^[42]. Note that TMS derivatives of alkene-triols, tetraols or pentaols obtained after NaBH4 reduction and derivatization of secondary oxidation products of di-, tri- or tetraunsaturated alkenones are too heavy and labile to be analyzed by GC-MS. The characterization of alkenone autoxidation products in sediments or phytodetritus with more adapted analytical techniques constitutes a very important challenge.

8. Carotenoids

Carotenoids, which are important antioxidant constituents of thylakoid membranes, play special roles in the protection of tissues against damage caused by light and oxygen ^[43]. These compounds can very efficiently quench ¹O2 by energy transfer (quenching), but also by chemical reaction (scavenging) ^[44]. They are also good scavengers of ROS ^[45]. The attack of β -carotene by ¹O2 affords β -carotene-5,8-endoperoxide (Scheme 10) ^[46]. If this compound is generally considered a useful early signal of ¹O2 production in plant leaves ^[47], it may be also formed during autoxidation of β -carotene ^[44] and is clearly not stable enough to serve as a viable environmental tracer. Unfortunately, the reaction of ¹O2 and ROS with carotenoids produces oxidation products that are not sufficiently stable and specific (production of similar compounds by enzymatic processes) ^[44] to be used as unequivocal indicators of type-II photosensitized oxidation or autoxidation of carotenoids in senescent phototrophic organisms and environmental samples.



Scheme 10. Reaction of ¹O2 with β -carotene.

Retrieved from https://encyclopedia.pub/entry/history/show/53164